How to Regenerate a PDF Report on Lubetrend



1. From the 'Results' Section, select the report(s) you wish to Regenerate:

2. Once the report(s) are selected, click on the 'Menu' button located at the bottom of the table and click on 'Generate PDF Reports'.

| Dashboard Results KPI | Recent Archived | | | | | | |
|--|---|--------|------|------------|------------|--------------|----------------|
| Pre-Reg Sample Labels Docs Help | Rec'd Date | - | From | 2022-12- | 08 ¢ To | 2023-12-09 | 4 Actions None |
| | □ 🗕 🖻 | • | / 🕮 | 足 🖿 | Lab No | Rec'd Date 🧅 | Flag |
| | | | | | x | | All 🗸 🗙 |
| | | | | 2 | 500601274 | 2023-11-06 | Т |
| | Export all to CSV Export selected to CSV Deselect all Exmail selected records to yourself | | | P h | 500601293 | 2023-11-06 | Т |
| | | | | E 🚹 | 500601266 | 2023-11-06 | Т |
| | | | | r 🖹 🐘 | 500601268 | 2023-11-06 | Т |
| | Generate P | amples | | 500601285 | 2023-11-06 | TI | |
| | Menu φ | | | | | | |

3. Lubetrend will ask you to confirm if you wish to continue, select 'Yes' to proceed.



4. Lubetrend will then Process the re-generated report(s) and confirm when complete.

You can then view them by selecting the Icon ²²¹ to download the report PDF file.